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This instruction for use HERMLE rotors supplement the operating manual for the centrifuge and do not replace it. Please also read the operating manual before starting up the rotors for the first time.


## Safety notes

### ⚠ ATTENTION: Fixing rotor, rotor lids and caps


- Please check before each run, if the rotor nut is tightened, rotor lids, caps and inserted carriers are correctly closed.
- If unusual noises occur when the centrifuge starts, the rotor, the rotor lid or cap may not be properly secured. Please stop centrifugation immediately by pressing the STOP key.

### ⚠ ATTENTION: Loading Rotors

- Rotors must be load symmetrically and with equal weight.

Picture 1

Picture 2

- The adapter may only be load with the appropriate vessels.
- Use always the same type of vessels and plates (material/density, volume weight). Always observe the max. g-force of the vessels and plates as indicated by the manufacturer
- Because of safety reason, every bucket with the same rating number has to be placed on the opposite site of the rotor. This provides for optimal operating characteristics and higher imbalance tolerance. The rating number is located on the top of each bucket.


Abbildung 3: Correct insertion of the buckets

Abbildung 4: Example rating number on the bucket

**ATTENTION: Aggressive Chemicals**

- Rotors are high-quality components that withstand extremely heavy loads. To ensure this stability it must be considered that aggressive chemicals can impair the stability of the rotor.
- Avoid the use of aggressive chemicals including strong weak alkali, strong acids, solutions with mercury, copper and other heavy metal ions, halogenated hydrocarbons, concentrated saline solutions and phenol.
- If the rotor, rotor lids and /or caps are contaminated with aggressive chemicals, immediately clean it using a neutral cleaning agent and then rinse it thoroughly with water. This applies to the rotor bores in particular.
- Rotor lids and caps made of PC and PP should be regularly checked of damages
- Replace rotor lids or caps which show cracks or milky stains immediately.

**Life time**

Rotors and rotor lids made of aluminium and stainless steel, have a operating time of **max. 7 years** from first use. Transparent rotor lids and caps made of PC and PP as well as rotors, tube racks and adapters of PP have a maximum operating time up to **3 years** from first use.

Condition for the operating time:

Proper use damage-free condition, recommended care.

- Please check before each run, if the accessories have damages. Replace any damaged accessories.
- Rotors, rotor lids, buckets, adapters or caps which showing any signs of corrosion or mechanical damage can't be use any longer.
- Do not use any accessories which are past their lifetime.
- Take care when inserting buckets and rotors, ensure that they do not become scratches.
- Protect the equipment from damage.

**Cleaning and disinfection**

For cleaning rotors and accessories (sleeves, adapters, buckets, rotor lids) use a mild cleaning agent. To prevent corrosion or material change, we recommend the cleaning at least once a week. Apart from this dirt must be removed immediately.

1. Remove the rotor from the centrifuge.
2. Use bottle brush to clean and disinfect the rotor bores.
3. Rotors and caps with sealing: Remove the seals and clean the groove
4. Rinse the rotors, rotor lid and adapter with clear water. Particular drilling of angle rotors.
5. For drying of the rotors and accessories set them on a towel. Place the angle rotors with bores down, to dry them to.
6. Dry the rotor cone with a soft, dry and lint-free cloth.
7. Do not grease the cone. Look for damage.
8. To prevent corrosion due to moisture between rotor and the motor shaft, the rotor should be at least once a week removed and the motor shaft should be cleaned with a clean, dry and lint-free cloth.

## Autoclaving

After Cleaning/disinfection the rotor can be autoclaved like described below.

The recommend time for autoclaving : 15-20 min at 121°C (1 bar)

- ⚠ ATTENTION: The sterilization time of 20 min. must not be exceeded. Sterilization again and again will cause reduction of the mechanical resistance of the plastic material. We recommend to replace the seals of aerosol-tight lids and caps after five autoclaving cycles.**

## Grease

1. After each cleaning and autoclaving process and if swing out rotors does not swing freely, apply a thin layer of Vaseline (order number 41-5215) to the rotor bolts (see picture 4).
2. Ensure that the rotor cross pivots and the bucket grooves are free from contamination.


Picture 3


Picture 4

3. After each run apply a thin layer of Vaseline (order number 41-5215) to the lid thread on fixed-angle rotors to prevent any damage (see picture 5 and 6). This prevents a twisting of the sealing during opening/closing. To ensure tightness, take care that the lid is tightly closed.

Picture 5

Picture 6

## Replace the seal

1. Use a blunt lever to lift the rubber seal from the nut. For this purpose use e.g. the round side of a paper clip (see picture 7). Take care and make sure that the seal will not damage.
2. Check the seal to make sure it is intact. Do not use any damaged, discolored or dirty seals.
3. Insert the seal in one location along the groove and use your fingers and spread it along the outer edge of the lid until the seal is fully seated in the groove.
4. Place the lid on the rotor and close the lid.

5. Remove the lid and check that the seal is seated correctly.

 If the seal is too long or too short, remove the seal from the groove and reinsert it.  
**ATTENTION: Faulty sealing when the gasket is handled incorrectly**  
**Uniformly insert the seal. Do not pull the seal lengthwise!**

Picture 7

### **Additional information for aerosol-tight rotors, lids and buckets**

The aerosol tightness of rotors, rotor lids, buckets and caps has been tested and certified by the „TÜV Nord CERT GmbH, Certification Body Consumer Products, Essen (Germany)“ in accordance with Annex AA of IEC 61010-2-020. The certificates can be downloaded on our webpage [www.hermle-labortechnik.de](http://www.hermle-labortechnik.de). Aerosol-tight rotors and buckets are marked with the label „aerosol-tight“.

 **ATTENTION: Autoclaving, mechanical stresses and contamination by chemicals or other aggressive solvents can impair the aerosol-tightness of the rotors and buckets.**

- Check the integrity of the seals of the aerosol-tight rotor lids or caps before each use.
- Use only aerosol-tight rotor lids or caps if the seals are undamaged and clean.
- Replace the seals of aerosol-tight lids and caps after five autoclaving cycles. .
- **Never** store aerosol-tight rotors or buckets closed.

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©HERMLE Labortechnik GmbH 2015
HERMLE Labortechnik GmbH